
LabJack

2023-April-03

1	Table of Contents	3
1.1	Overview	3
1.2	Supported models	3
1.2.1	T4	3
1.2.2	T7	5
1.2.3	T7-PRO	6
1.2.4	T8	7
1.3	Installation	8
1.4	Driver architecture	8
1.5	Configuration	8
1.6	Top-level OPI screen	9
1.7	Databases	9
1.7.1	Device Functions	9
1.7.2	Analog Input Functions	11
1.7.3	Analog Output Functions	11
1.7.4	Digital I/O Functions	13
1.7.5	Waveform Digitizer Functions	14
1.7.6	Waveform Generator Functions	15

author Mark Rivers, University of Chicago

An [EPICS](#) driver for the T-series devices from [LabJackCorp](#). These multi-function devices support support analog input, thermocouple input, analog output, and binary I/O.

1.1 Overview

This is an [EPICS](#) driver for the T-series devices from [LabJackCorp](#). These multi-function devices support support analog input, thermocouple input, analog output, and binary I/O.

The driver is written in C++, and consists of a class that inherits from [asynPortDriver](#), which is part of the EPICS [asyn](#) module.

The driver is written to be general, so that it can be used with any LabJack T-series module. The T-series devices all use Modbus for the low-level communications, and use a consistent set of Modbus register addresses, so the code is largely model-independent. The driver does require small modifications to be be used with a new model.

1.2 Supported models

All LabJack T-series modules have both Ethernet and USB communications. The following models are currently supported.

1.2.1 T4

This module costs \$245 and has the following features:

- 4 single-ended analog inputs
 - +- 10V range
 - 12-bit
 - Up to 50 kHz total streaming input rate, i.e. 1 channel at 50 kHz, 2 channels at 25 kHz, etc.
- Up to 8 additional analog inputs (0-2.5 V range, 12-bit). These can also be configured as digital I/O bits.
- 2 analog outputs



Fig. 1: Photo of LabJack-T4

- 0-5V range
- 10-bit
- Up to 50 kHz streaming output rate
- 20 digital I/O bits
 - Each configurable as input or output.
 - 8 of these can be configured as 0-2.5V analog inputs.

More information can be found in the [T4 product description](#).

1.2.2 T7



Fig. 2: Photo of LabJack-T7

This module costs \$520 and has the following features:

- 14 single-ended / 7 differential analog inputs
 - Programmable range (+-10V, +-1V, +-0.1V, +-0.01V)
 - 16-bit
 - Up to 100 kHz total streaming input rate, i.e. 1 channel at 100 kHz, 2 channels at 50 kHz, etc.

- 2 analog outputs
 - 0-5V range
 - 12-bit
 - Up to 100 kHz streaming output rate
- 23 digital I/O bits - Each configurable as input or output.

More information can be found in the [T7 product description](#).

1.2.3 T7-PRO




Fig. 3: Photo of LabJack T7-PRO

This module costs \$750 and is the same as the LabJack T7 with the following additional features:

- Each analog input can be software switched between 16-bit and 24-bit ADC, trading off resolution for speed.
- In 24-bit mode the ADCs support 9 types of thermocouple inputs.
- WiFi communications (in addition to standard Ethernet and USB).

More information can be found in the [LabJack T7-PRO product description](#).

1.2.4 T8




Fig. 4: Photo of LabJack T8

This module costs \$1,400 and has the following features:

- 8 differential analog inputs
 - 11 input ranges from $\pm 11\text{V}$ to $\pm 0.15\text{V}$
 - 24-bit sigma/delta ADC
 - Supports 9 types of thermocouple inputs
 - Up to 40 kHz streaming input per channel, i.e. up to 320 kHz total scanning rate.
- 2 analog outputs
 - 0-10V range
 - 16-bit
 - Up to 40 kHz streaming output rate
- 23 digital I/O bits, each configurable as input or output.

More information can be found in the [LabJack T8 product description](#).

1.3 Installation

The EPICS LabJack module uses the [LJM library](#) from LabJack. It runs on Linux and Windows.

Most Linux versions should be supported. It has been tested on Centos 7. The EPICS module includes the LJM header and library files, so ideally LJM would not need to be installed locally on Linux. However, LJM uses configuration files which it installs in `/usr/local/share/LabJack/LJM`, so the LJM package does need to be installed, and this requires root privilege. I have asked LabJack support if it is possible for those files to be located in a directory that does not require root permission to write to.

On Windows the LJM library package needs to be installed to run the IOC.

1.4 Driver architecture

The driver has a polling thread that periodically reads the state of the digital I/O bits and the values of the analog inputs. If the waveform digitizer or waveform generator functions are active it polls the status of those as well. The delay time at the end of the polling cycle can be controlled via an EPICS PV. The actual poll cycle time, including the delay, is reported in an EPICS PV.

The digital I/O are normally set to `SCAN=I/O Intr` so that they change state quickly when the poller reads them.

The analog inputs can be set to either periodic or I/O Intr scan. I/O Intr scan allows more frequent updates at the expense of CPU load and Channel Access traffic.

When reading analog inputs in thermocouple mode the inputs are actually read in volts, and the conversion to temperature is done in software. This uses the cold junction temperature read from the device, and the temperature conversion function in the LJM library. This allows temperature inputs to be scanned with the waveform digitizer function, which is not possible if the temperature conversions are performed on the device itself.

1.5 Configuration

The following lines are needed in the EPICS startup script for the LabJack driver.

```
## Configure port driver
# LabJackConfig(portName,          # The name to give to this asyn port driver
#               uniqueID,         # The IP name, IP address, or serial number of the
↳LabJack module.
#               maxInputPoints,   # Maximum number of input points for waveform
↳digitizer
#               maxOutputPoints) # Maximum number of output points for waveform
↳generator
LabJackConfig("LJT7_1", "gse-labjack1.cars.aps.anl.gov", 2048, 2048)
```

The `uniqueID` is a string that identifies the device to be controlled. It can be any of the following:

- A fully qualified domain name with periods, e.g. `gse-labjack1.cars.aps.anl.gov`. The periods are needed to distinguish an IP name from a serial number.
- An IP address, e.g. `10.54.160.72`.
- A module serial number, e.g. `470029169`.

The LabJack module comes with example `iocBoot/` directories that contain example startup scripts and example substitutions files for each model.

1.6 Top-level OPI screen

The following is the top-level medm screen for the LabJack T7 and T7-PRO modules.

The screenshot shows the LabJack T7 top-level OPI screen with the following sections:

- Header:** LabJack LabJackT7_1:
- System Information:**
 - Model name: T7
 - Serial number: 470029169
 - LJM version: 1.220000
 - Last error: No error
 - Firmware version: 1.029200
 - Temperature (K): 300.369
 - Driver version: 3.0.0
 - Poll sleep time (ms): 10.0
 - Poll cycle time (ms): 15.2
 - Reset device on exit: **Reset!**
- Digital I/O:**
 - FIO 0-7:** Value 0xf. Channels 0-3 are red (In), 4-7 are green (In).
 - EIO 0-7:** Value 0xff. Channels 8-15 are red (In).
 - CIO 0-3:** Value 0xf. Channels 16-19 are red (In), 20-22 are green (In).
 - MIO 0-2:** Value 0x7. Channels 20-22 are green (In).
- Analog Input:**
 - Channels 0-13 with values: 1.1352, 1.1350, 0.0547, 0.0379, -9999.9999, -9999.9999, -9999.9999, -9999.9999, -9999.9999, -9999.9999, -9999.9999, -9999.9999, -9999.9999, -9999.9999.
 - Buttons: Strip charts, Configure, Waveform digitizer.
- Analog Output:**
 - Channels 0-1 with values: 1.1347, 1.1346.
 - Buttons: Waveform generator, Configure, LJTICK-DACS.

Fig. 5: LabJack_T7.adl

While this screen is nominally specific to the T7 and T7-PRO, it can be used for any model. On the T8 analog inputs 8-13 do not apply. On the T4 analog inputs 12-13 do not apply, nor do digital I/O bits 20-22.

1.7 Databases

The following tables list the database template files that are used with the LabJack modules.

1.7.1 Device Functions

These are the records defined in LabJack_device.template. These records provide device-wide information and control. This database is loaded once for each LabJack device.

EPICS record name	EPICS record type	asyn interface	drvInfo string	Description
\$(P)ModelName	mbbi	asynInt32	MODEL_NAME	The model name. mbbi values and strings are <ul style="list-style-type: none"> • 0="T4" • 1="T7" • 2="T7-Pro" • 3="T8".
\$(P)FirmwareVersion	asynVersion	asynOctet	FIRMWARE_VERSION	Firmware version.
\$(P)SerialNumber	asynInt32	asynOctet	SERIAL_NUMBER	Serial number.
\$(P)DeviceTemperature	asynFloat64	asynFloat64	DEVICE_TEMPERATURE	Temperature. This is used as the cold junction reference temperature for thermocouple measurements. It has SCAN="5 second" which is fast enough for this slowly varying value.
\$(P)LJMVersion	asynVersion	asynOctet	LJM_VERSION	Version of the LabJack LJM library.
\$(P)DriverVersion	asynVersion	asynOctet	DRIVER_VERSION	Version of the EPICS driver.
\$(P)LastErrorMessage	asynOctet	asynOctet	LAST_ERROR_MESSAGE	Message from the driver. This includes a timestamp.
\$(P)PollSleepMS	asynFloat64	asynFloat64	POLL_SLEEP_MS	The number of milliseconds to sleep at the end of each poll cycle.
\$(P)PollTimeMS	asynFloat64	asynFloat64	POLL_TIME_MS	The actual number of milliseconds to execute the poll cycle, including the sleep.
\$(P)AiAllScaling	asynFloat64	asynFloat64	ANALOG_IN_SCALING	Set the scaling for all analog input channels.
\$(P)AiAllResolution	asynInt32	asynInt32	ANALOG_IN_RESOLUTION	Resolution. Lower resolution result in lower noise and longer ADC conversion time. Resolution 0 is the default resolution for that model. <ul style="list-style-type: none"> • The T4 supports resolutions 1-5. • The T7 supports resolutions 1-8. • The T7-PRO supports resolutions 1-12. 1-8 use the 16-bit ADC and 9-12 use the 24-bit ADC. When running the waveform generator on the T7-PRO this must be set to values between 1-8, i.e. 16-bit ADC. The driver will set this automatically when starting the waveform generator if it is outside the allowed range. • The T8 supports resolutions 1-16. However, it is recommended to use the default resolution and change the SamplingRate to control the resolution vs speed tradeoff.
\$(P)AiSamplingRate	asynFloat64	asynFloat64	ANALOG_IN_SAMPLING_RATE	Rate of the ADC in Hz. <ul style="list-style-type: none"> • It applies to the T8 only. • Recommended range is 100 to 10000 Hz. • Lower rates do more filtering in the ADC, reducing noise at the expense of speed. • Increasing the sampling rate will increase the noise in each reading. • However, since the analog input records use the devAsynFloat64Average device support, increasing the rate can increase the number of samples averaged in the EPICS device support in a fixed period of time, provided it is not limited by PollSleepMS. • Because of this averaging in device support, increasing the sampling time from 100 Hz to 1000 Hz can actually result in a small decrease in noise. • The maximum rate that the values can be read from the device with PollSleepMS=0 is about 2000/s, so increasing the SamplingRate beyond 2000 will not result more averaging in EPICS device support.

1.7.2 Analog Input Functions

These are the records defined in LabJack_ai.template. This database is loaded once for each analog input channel

EPICS record name	EPICS record type	asyn interface	drvInfo string	Description
\$(P)Ai\$(N)ai	asynFloat64	ANALOG	ANALOG	ANALOG Input value. This is polled in the driver, so either period or I/O Intr scanning can be used.
\$(P)AiEnable\$(N)	asynInt32	ANALOG	ANALOG	ENABLE For this analog input channel. Disabled inputs are not read by the poller. Unconnected inputs should be disabled to improve accuracy on active channels and to reduce the polling time.
\$(P)AiModel\$(N)	asynInt32	ANALOG	ANALOG	INPUT MODE for this analog input channel. Choices are Volts and 9 different thermocouple types.
\$(P)AiTempUnits\$(N)	asynInt32	TEMPERATURE	TEMPERATURE	TEMPERATURE Units for this analog input channel if a thermocouple mode is selected. Choices are “K”, “C”, and “F”.
\$(P)AiDiff\$(N)	asynInt32	ANALOG	ANALOG	MODE Single-Ended” or “Differential” input mode on the T7 and T7-PRO. The T4 is always single-ended and the T8 is always differential. The driver constructs the strings and values based on the model.
\$(P)AiRange\$(N)	asynInt32	ANALOG	ANALOG	SELECT RANGE input range for this analog input channel. <ul style="list-style-type: none"> On the T4 the range is fixed at +-10V on channels 0-3 and 0-2.5 on channels 4-11. On the T7 the range choices are +-10V, +-1V, +-0.1V, and +-0.01V. On the T8 there are 11 ranges from +-11V to +-0.15V. The driver constructs the strings and values based on the model.
\$(P)AiResolution\$(N)	asynInt32	ANALOG	ANALOG	SELECT RESOLUTION resolution for this analog input channel. High values of resolution result in lower noise and longer ADC conversion time. <ul style="list-style-type: none"> Resolution 0 is the default resolution for that model. The T4 supports resolutions 1-5. The T7 supports resolutions 1-8. The T7-PRO supports resolutions 1-12. 1-8 use the 16-bit ADC and 9-12 use the 24-bit ADC The T8 supports resolutions 1-16. However, these are automatically selected by the Range, and this record has no effect?

The following is the medm screen for controlling the analog input configuration records.

While this screen is nominally specific to the T7 and T7-PRO, it can be used for any model. On the T8 analog inputs 8-13 do not apply and the inputs are always in Differential mode. On the T4 analog inputs 12-13 do not apply. Thermocouples only work well with the T7-PRO in 24-bit mode (resolutions 9-12), or with the T8 with low-voltage ranges. They do not work well with the T4 or T7.

1.7.3 Analog Output Functions

These are the records defined in LabJack_Ao.template. This database is loaded once for each analog output channel

LabJack_T7_AiSetup.adl@corvette

Analog Input Configuration LabJackT8_1:

Chan.	Enable	Mode	Range	Resol.	T	Low	High	Prec	SE/Diff.
0	Disable	Volts	± 0.15V	Default	C	-0.0003	-0.0002	8	Differential
1	Enable	Volts	± 1.2V	Default	C	-0.0001	-0.0001	6	
2	Disable	Volts	± 11.0V	Default	K	0.0000	0.0000	6	Differential
3	Disable	Volts	± 11.0V	Default	C	0.0000	0.0000	4	
4	Disable	Volts	± 11.0V	Default	K	0.0000	0.0000	4	Differential
5	Disable	Volts	± 11.0V	Default	K	0.0000	0.0000	4	
6	Disable	Volts	± 11.0V	Default	K	0.0000	0.0000	4	Differential
7	Disable	Volts	± 0.18V	Default	K	0.0000	0.0010	6	
8	Disable	Volts	N.A.	N.A.	K	0.0000	0.0000	4	N.A.
9	Disable	Volts	N.A.	N.A.	K	0.0000	0.0000	4	
10	Disable	Volts	N.A.	N.A.	K	0.0000	0.0000	4	N.A.
11	Disable	Volts	N.A.	N.A.	K	0.0000	0.0000	4	
12	Disable	Volts	N.A.	N.A.	K	0.0000	0.0000	4	N.A.
13	Disable	Volts	N.A.	N.A.	K	0.0000	0.0000	4	

Settling time (us) (0=auto)
 Resolution (all channels)
 Sampling rate (Hz) (T8 only)

NOTE: on T7 and T7-Pro in differential mode, channels are paired 0+ 1-. 2+ 3-. etc.

Fig. 6: LabJack_T7_AiSetup.adl

EPICS record name	EPICS record type	asyn interface	drvInfo string	Description
\$(P)\$R)ao	ao	asynFloat64	ANALOG	Analog output value.
\$(P)\$R)TweakVal	TweakVal	N.A.	N.A.	The amount by which to tweak the out when the Tweak record is processed.
\$(P)\$R)TweakUp	TweakUp	N.A.	N.A.	Tweaks the output up by TweakVal.
\$(P)\$R)TweakDown	TweakDown	N.A.	N.A.	Tweaks the output down by TweakVal.

The following is the medm screen for configuring the analog output records. The drive limits can be more restrictive than the full output range of the analog outputs.

Channels 0 and 1 are the on-board DACs on all models. The range is 0-5V for the on-board DACs on the T4, T7, and T7-PRO, and 0-10V on the T8.

Channels 2 and above are the LJTick DACs that can be optionally installed on any model. The range on the LJTick DACs is +-10V.




Fig. 7: LabJack_T7_AoSetup.adl

The following is the medm screen for controlling the optional LJTick DACs.




Fig. 8: LabJack_LJTick_DAC2.adl

1.7.4 Digital I/O Functions

These are the records defined in LabJack_binary.template and LabJack_biWord.template.

EPICS record name	EPICS record type	asyn interface	drvInfo string	Description
\$(P)Bi\$(N)bi		asynUInt32DIGITAL	DIGITAL	Digital input value. The MASK parameter in the INP link defines which bit is used. The binary inputs are polled by the driver poller thread, so these records should have SCAN="I/O Intr".
\$(P)Bo\$(N)bo		asynUInt32DIGITAL	DIGITAL	Digital output value. The ADDR parameter in the INP link defines which bit is used.
\$(P)Bd\$(N)bo		asynUInt32DIGITAL	DIGITAL	DIRECTION of this I/O line, "In" (0) or "Out" (1). The MASK parameter in the INP link defines which bit is used.
\$(P)\$R	longin	asynUInt32DIGITAL	DIGITAL	Digital input value as a word, rather than individual bits. The ADDR parameter in the INP link defines which word is read. 0=DIO (bits 0-23), 1=FIO (bits 0-7), 2=EIO (bits 8-15), 3=CIO (bits 16-19), and 4=MIO (bits 20-22). The binary inputs are polled by the driver poller thread, so these records should have SCAN="I/O Intr".

1.7.5 Waveform Digitizer Functions

These records are defined in the following files: - LabJack_waveformDig.template. This database is loaded once per module. - LabJack_waveformDigN.template. This database is loaded for each digitizer input channel.

EPICS record name	EPICS record type	asyn interface	drvInfo string	Description
\$(P)WaveDigNumPoints	asynInt32	32	WAVEDIG_NUMPOINTS	Number of points to digitize. This cannot be more than the value of maxInputPoints that was specified in LabJackConfig.
\$(P)WaveDigFirstChan	asynInt32	32	WAVEDIG_FIRST_CHAN	First channel to digitize, 0-13.
\$(P)WaveDigNumChans	asynInt32	32	WAVEDIG_NUMCHANS	Number of channels to digitize. 1-14. The maximum valid number is 13-FirstChan+1.
\$(P)WaveDigWaveWF	asynFloat32	64	WAVEDIG_TIMEASWF	Time as waveform. These values are calculated when Dwell or NumPoints are changed. It is typically used as the X-axis in plots.
\$(P)WaveDigCurrentPoints	asynInt32	32	WAVEDIG_CURRENT_POINTS	The current points being collected. This does not always increment by 1 because the device can transfer data in blocks.
\$(P)WaveDigDwell	asynFloat64	64	WAVEDIG_DWELL	Dwell time per point in seconds. The minimum time depends on the device type and NumChans.
\$(P)WaveDigDwellActual	asynFloat64	64	WAVEDIG_DWELL_ACTUAL	Actual dwell time per point in seconds. This may differ from the requested Dwell because of clock granularity in the device.
\$(P)WaveDigTotalTime	asynFloat64	64	WAVEDIG_TOTAL_TIME	Total time to digitize NumChans*NumPoints.
\$(P)WaveDigResolution	asynInt32	32	WAVEDIG_RESOLUTION	Resolution to use for all channels during the scan. The choices are model-dependent and are set by the driver.
\$(P)WaveDigSettlingTime	asynFloat64	64	WAVEDIG_SETTLING_TIME	Settling time in microseconds to use for all channels during the scan. 0 selects the device default.
\$(P)WaveDigExtTrigger	asynInt32	32	WAVEDIG_EXT_TRIGGER	External trigger, “Internal” (0) or “External” (1). NOTE: NOT YET IMPLEMENTED.
\$(P)WaveDigExtClock	asynInt32	32	WAVEDIG_EXT_CLOCK	External clock, “Internal” (0) or “External” (1). If External is used then the Dwell record does not control the digitization rate, it is controlled by the external clock. However Dwell should be set to approximately the correct value if possible, because that builds the time axis for plotting. NOTE: NOT YET IMPLEMENTED.
\$(P)WaveDigAutoRestart	asynInt32	32	WAVEDIG_AUTO_RESTART	“Auto-Restart” (0) and “Enable” (1). This controls whether the driver automatically starts another acquire when the previous one completes.
\$(P)WaveDigRun	asynInt32	32	WAVEDIG_RUN	“Stop” (0) and “Run” (1). This starts and stops the waveform digitizer. It will automatically stop when the requested number of samples have been acquired.
\$(P)WaveDigReadWF	asynInt32	32	WAVEDIG_READ_WF_DONE	“Read Done” (0) and “Read” (1). This reads the waveform data from the device buffers into the waveform records. Note that the driver always reads device when acquisition stops, so for quick acquisitions this record can be Passive. To see partial data during long acquisitions this record can be periodically processed.
\$(P)VoltWFs(N)	asynFloat64	64	WAVEDIG_VOLT_WF	Waveform record contains the digitizer waveform data for channel N. This record has scan=I/O Intr, and it will process whenever acquisition completes, or whenever the ReadWF record above processes. The data are in volts or temperature units.

This is the medm screen that controls the waveform digitizer.

This is a plot of the waveform digitizer waveform for analog input 0 capturing 1000 samples/s (Dwell=0.001).

1.7.6 Waveform Generator Functions

These records are defined in the following files: - LabJack_waveformGen.template. This database is loaded once per module. - LabJack_waveformGenN.template. This database is loaded for each waveform generator output channel.




Fig. 9: LabJackWaveDig.adl




Fig. 10: LabJackWaveDigPlot.adl

EPICS record name	EPICS record type	asyn interface	drvInfo string	Description
\$(P)WaveGenNumPoints	Int32	asynInt32	WAVEGEN	Number of points output waveform. The value of this record is equal to UserNumPoints if user-defined waveforms are selected, or IntNumPoints if internal predefined waveforms are selected.
\$(P)WaveGenUserNumPoints	Int32	asynInt32	WAVEGEN	Number of user-defined output waveforms. This cannot be more than the value of maxOutputPoints that was specified in LabJackConfig.
\$(P)WaveGenIntNumPoints	Int32	asynInt32	WAVEGEN	Number of internal predefined output waveforms. This cannot be more than the value of maxOutputPoints that was specified in LabJackConfig.
\$(P)WaveGenUserDwell	Float32	asynFloat32	WAVEGEN	User-based dwell time for user-defined waveforms. These values are calculated when UserDwell or UserNumPoints are changed. It is typically used as the X-axis in plots.
\$(P)WaveGenIntDwell	Float32	asynFloat32	WAVEGEN	Internal-based dwell time for internal predefined waveforms. These values are calculated when IntDwell or IntNumPoints are changed. It is typically used as the X-axis in plots.
\$(P)WaveGenCurrentPoints	Int32	asynInt32	WAVEGEN	Current number of points being output. This does not always increment by 1 because the device can transfer data in blocks.
\$(P)WaveGenFrequency	Float64	asynFloat64	WAVEGEN	Output frequency (waveforms/second). The value of this record is equal to UserFrequency if user-defined waveforms are selected, or IntFrequency if internal predefined waveforms are selected.
\$(P)WaveGenDwell	Float64	asynFloat64	WAVEGEN	Output dwell time or period (seconds/sample). The value of this record is equal to UserDwell if user-defined waveforms are selected, or IntDwell if internal predefined waveforms are selected.
\$(P)WaveGenDwellActual	Float64	asynFloat64	WAVEGEN	The actual dwell time. This can be different from the requested dwell time (WaveGenDwell) because of the granularity of the device clock.
\$(P)WaveGenUserDwell	Float64	asynFloat64	WAVEGEN	The user dwell time or period (seconds/sample) for user-defined waveforms. This record is automatically changed if UserFrequency is modified.
\$(P)WaveGenIntDwell	Float64	asynFloat64	WAVEGEN	The internal dwell time or period (seconds/sample) for internal predefined waveforms. This record is automatically changed if IntFrequency is modified.
\$(P)WaveGenUserFrequency	N/A	N/A	N/A	The output frequency (waveforms/second) for user-defined waveforms. This record computes UserDwell and writes to that record. This record is automatically changed if UserDwell is modified.
\$(P)WaveGenIntFrequency	N/A	N/A	N/A	The output frequency (waveforms/second) for internal predefined waveforms. This record computes IntDwell and writes to that record. This record is automatically changed if IntDwell is modified.
\$(P)WaveGenTotalTime	Float64	asynFloat64	WAVEGEN	The total time to output the waveforms. This is WaveGenDwellActual*NumPoints.
\$(P)WaveGenExtTrigger	Int32	asynInt32	WAVEGEN	TEXT_TRIGGER, "Internal" (0) or "External" (1). NOTE: NOT YET IMPLEMENTED,
\$(P)WaveGenExtClock	Int32	asynInt32	WAVEGEN	TEXT_CLOCK, "Internal" (0) or "External" (1). If External is used then the Dwell record does not control the output rate, it is controlled by the external clock. However Dwell should be set to approximately the correct value if possible, because that controls the time axis on the plots. NOTE: NOT YET IMPLEMENTED.
\$(P)WaveGenContinuous	Int32	asynInt32	WAVEGEN	WAVEGEN_CONTINUOUS, "One-shot" (0) or "Continuous" (1). This controls whether the device stops when the output waveform is complete, or immediately begins again at the start of the waveform.
\$(P)WaveGenRun	Int32	asynInt32	WAVEGEN	WAVEGEN_RUN are "Stop" (0) and "Run" (1). This starts and stops the waveform generator. In one-shot mode the waveform generator stops automatically when all of the samples have been output.
\$(P)WaveGenUserWF\$(N)	Float64 Array	asynFloat64Array	WAVEGEN	This user waveform record contains the user-defined waveform generator data for channel N. The data are in volts. These data are typically

This is the medm screen for the waveform generator.




Fig. 11: LabJackWaveGen2.adl

This is a plot of an internal predefined sin wave waveform.




Fig. 12: Plot of an internal predefined waveform (sin wave)